

Air Shower



Accredited for compliance with ISO/IEC 17025-testing



MODELS	AS-1P1S	AS-1P2S	AS-2P2S
External Size (W*D*H)	1300 x 1000 x 2050 mm	1500 x 1000 x 2050 mm	1500 x 2000 x 2050 mm
Shower Area Size (W*D*H)	800 x 890 x 1930 mm	800 x 890 x 1930 mm	800 x 890 x 1930 mm
HEPA Filter	99.999% efficiency at 0.3µm		
Spraying Wind Speed	25m/s		
Air Shower Time	0~99s (adjustable)		
No. of SS Nozzles	6	12	24
Fluorescent Lamp	12W*1	12W*1	12W*2
Cabinet Material (Optional)	1. Full cold-rolled steel with powder coating 2. Internal stainless steel, external cold-rolled steel 3. Full stainless steel		
Door Material (Optional)	1. Full cold-rolled steel with powder coating 2. Stainless steel		
Chassis	Stainless steel		
Power Supply	AC220V±10%, 50/60Hz 110V±10%, 60Hz		
Consumption	1000W	1500W	2000W
Gross Weight	300kg	400kg	700kg
Package size (W*D*H)	1400 x 1150 x 2200mm	1660 x 1150 x 2200mm	1600 x 1650 x 2200mm



OPIRA BRISBANE
32 DIVIDEND STREET
MANSFIELD, QLD 4122

OPIRA MELBOURNE
25 GRAHAM ROAD
CLAYTON SOUTH, VIC 3169

OPIRA PERTH
9 EARLSTON PLACE
BOORAGOON, WA 6154

OPIRA SYDNEY
20 DUKE ST
FORESTVILLE, NSW 2087

CONTACT US

1300 157 969

SOLUTIONS@OPIRA.COM.AU