

Dispensing Booth



Accredited for compliance with ISO/IEC 17025-testing

A dispensing booth is a type of partial purifying equipment for filling, refilling, weighing and sampling raw materials and compounds. Supplied with a HEPA filter, which prevents airborne specks of dusts by the down draught technique.



Our range of dispensing booths features:

- Designed and manufactured to comply with AS 2252.6-2011 Controlled environments Clean workstations - Design, install and use
- Differential pressure gauge - Dwyer make Magnehelic gauges that measure the pressure drop HEPA filter 0-500 Pa
- LCD Display, microprocessor control system
- Soft wall material - PVC anti-static dustproof curtain



OPIRA BRISBANE
32 DIVIDEND STREET
MANSFIELD, QLD 4122

OPIRA MELBOURNE
25 GRAHAM ROAD
CLAYTON SOUTH, VIC 3169

OPIRA PERTH
9 EARLSTON PLACE
BOORAGOON, WA 6154

OPIRA SYDNEY
20 DUKE ST
FORESTVILLE, NSW 2087

CONTACT US

1300 157 969

SOLUTIONS@OPIRA.COM.AU

Technical Specifications

MODELS	OPIRA-DB-1200	OPIRA-DB-1800	OPIRA-DB-2400
Clean Level	ISO 5 (Class 100), Class A		
HEPA Filter	99.999% efficiency at 0.3µm		
Air Velocity	Supply air volume ≤7500m ³ /h, Exhaust air volume ≤2250m ³ /h; adjustable		
Pressure Gauge	3 pcs		
UV Lamp - emission of 253.7 nanometers, with timer	30W*1	40W*1	30W*2
LED Lamp ≥500lux	14W*4	16W*4	16W*4 & 14W*2
Power Supply	AC220V±10%, 50/60Hz 110V±10%, 60Hz		
Consumption	600W	1000W	1700W
Sockets	2 pcs		
Work Area Size (W*D*H)	1200 x 1200 x 2000mm	1800 x 1800 x 2000mm	2400 x 1800 x 2000mm
External Size (W*D*H)	1300 x 1700 x 2400mm	1900 x 2300 x 2400mm	2500 x 2400 x 2400mm
Gross Weight	420kg	750kg	840kg
Package Size (W*D*H)	1440 x 2540 x 1980mm	2020 x 760 x 2420mm 2020 x 660 x 495mm 1950 x 530 x 2160mm	2640 x 960 x 2380mm 2640 x 780 x 680mm 2640 x 560 x 2250mm

Customised design for clean level, size or material is available.



Our products have been designed and tested to comply with all applicable ACMA regulatory requirements and bear the RCM mark.



OPIRA BRISBANE
32 DIVIDEND STREET
MANSFIELD, QLD 4122

OPIRA MELBOURNE
25 GRAHAM ROAD
CLAYTON SOUTH, VIC 3169

OPIRA PERTH
9 EARLSTON PLACE
BOORAGOON, WA 6154

OPIRA SYDNEY
20 DUKE ST
FORESTVILLE, NSW 2087

CONTACT US

1300 157 969

SOLUTIONS@OPIRA.COM.AU