

Passbox



Accredited for compliance with ISO/IEC 17025-testing



Our products have been designed and tested to comply with all applicable ACMA regulatory requirements and bear the RCM mark.

TECHNICAL SPECIFICATIONS

Model	ASPB-01	ASPB-01	ASPB-01	PB-01	PB-02	PB-03
External size (W*D*H) mm	660*570*1000	760*670*1100	860*770*1200	660*570*580	760*670*680	860*770*780
Internal size (W*D*H) mm	500*500*500	600*600*600	700*700*700	500*500*500	600*600*600	700*700*700
Material	304 stainless steel			304 stainless steel		
Interlock	Electronic			Electronic		
UV lamp	30W*1			30W*1		
HEPA filter	Efficiency≥ 99.999% at 0.3µm			Not applicable		
Clean air speed	0.36~0.54m/s			Not applicable		
Gross weight	75kg	100kg	130kg	35kg	40kg	50kg
Package size	700*770*1200	800*870*1300	900*970*1400	700*700*780	800*800*880	900*900*980



OPIRA BRISBANE
32 DIVIDEND STREET
MANSFIELD, QLD 4122

OPIRA MELBOURNE
25 GRAHAM ROAD
CLAYTON SOUTH, VIC 3169

OPIRA PERTH
9 EARLSTON PLACE
BOORAGOON, WA 6154

OPIRA SYDNEY
20 DUKE ST
FORESTVILLE, NSW 2087

CONTACT US

1300 157 969

SOLUTIONS@OPIRA.COM.AU