


Protecting First Responders & Patients against Airborne Pathogens

VP900 INTERCEPTOR
Air Disinfection Unit for Ambulance

Portable unit that uses germicidal UV to disinfect the air. Designed to be mounted in the interior of vehicle (under a seat or on a wall).

The VP900 intercepts airborne pathogens and contaminants to insure drivers and patients safety.



BENEFITS ▾

- Destroys airborne pathogens
- Reduces the risk of contamination during transportation
- Easy installation

FEATURES ▾


- 12V power supply
- Designed to be incorporated into emergency vehicles
- Shock-resistant construction made of steel and aluminum
- Unit dimensions: 43 cm x 28 cm x 13cm

WARRANTY ▾

- Lamp: 2 years
- Ballast: 3 years



UV 'J' Lamp



Flexible installations



Eliminates odors and chemical contaminants



Eliminates germs



Eliminates volatile organic compounds (VOCs)